REF IS-6010 	IDS-iSYS CARTRIDGE CHECK SYSTEM (CCS)	 immundiagnostic systems
Package Insert	IN VITRO DIAGNOSTIC	CE

## Intended use

### For In Vitro Diagnostic Use

The IDS-iSYS Cartridge Check System (CCS) is a quality control cartridge used to verify the performance of various modules used for performing immunoassays on the IDS-iSYS Multi-Discipline Automated Analyser.

## Summary and Explanation

The CCS cartridge contains various reagents dedicated to the verification of the pipetting, incubation, washing and reading modules involved in performing immunoassays on the IDS-iSYS analyser.

## Method Description

The analytical principles of the immunoassays employed on the IDS-iSYS analyser are based on chemiluminescent (ICMA, LIA) signals and use antibodies or antigens, either immobilized or labelled by biotin, or labelled by an acridinium ester derivative (DMAE). The antibodies or antigens are used in two-site (sandwich) or one-site (competitive) immunoassays to determine the concentrations of various analytes.

After a washing step, the light emitted by the DMAE label bound to magnetic particles is directly or inversely proportional to the amount of analyte present in the calibrators or the samples.

The testing protocols ran on the IDS-iSYS analyser to assess its performance level are as follows:

**CCS1**: Test protocol to check the performance of the distribution and luminometer modules. High RLU value expected.

**CCS2**: Test protocol to check the performances of the distribution and luminometer modules. Low RLU expected.

**CCS3**: Test protocol to check the signal background after washing.

**CCS4**: Test protocol to check the signal background for use in troubleshooting.

**CCSB**: Test protocol for use in troubleshooting, to be used under direction of Field Service Support.

## Warnings and Precautions

The IDS-iSYS Cartridge Check System (CCS) is for in vitro diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Package Insert. IDS Limited will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of non-compliance with the instructions provided.

**CAUTION:** This kit contains material of human and animal origin. Handle kit reagents as if capable of transmitting an infectious agent.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

### Human material

Human material used in the preparation of this product has been tested by FDA recommended assays for the presence of antibody to Human Immunodeficiency Virus (HIV I and II), Hepatitis B surface antigen, antibody to

Hepatitis C, and found negative. As no test can offer complete assurance that infectious agents are absent, the reagents should be handled in accordance at Biosafety Level 2.

### Sodium Azide

Some reagents in this kit contain sodium azide <0.1% (w/w) which may react with lead, copper or brass plumbing to form highly explosive metal azides. On disposal, flush with large volumes of water to prevent azide build up.

## Handling precautions

The reagents provided in the kit are ready to use. When a new cartridge is loaded onto the analyser, the magnetic particle container requires mixing by brisk rotation motion. This will resuspend the magnetic particles that have settled during shipment. It is very important to avoid any foam formation.

## Shelf Life and Storage of Reagents

Store the cartridges in a vertical position at 2 to 8 °C. Do not freeze the cartridge.

Storage temperature	Reagent shelf life
Before opening at 2 – 8 °C	Refer to expiration date
After opening at 2 – 8 °C	9 weeks
On board the IDS-iSYS analyser	8 weeks

## Procedure

### Materials Provided

#### CCS Cartridge

A CCS kit is comprised of 1 cartridge sufficient to run the tests presented for a 2 month period:

REAG	1	MP	2 ml
REAG	2	CONJ	20 ml
REAG	3	Ab-BIOT	9 ml
REAG	4	DIL	9 ml

**REAG 1** – Magnetic particles coated with streptavidin in phosphate buffer containing BSA and sodium azide as preservative (<0.1%).

**REAG 2** : Anti-hIgG mouse monoclonal antibody labeled with an acridinium ester derivative in buffer containing a stabilizing agent and sodium azide as preservative (<0.1%).

**REAG 3** : Human immunoglobulin conjugated to biotin in a buffer containing a stabilizing agent and sodium azide as preservative (<0.1%).

**REAG 4** : Phosphate buffer containing detergent and sodium azide as preservative (<0.1%).

### Materials Required But Not Provided

#### IDS-iSYS Multi-Discipline Automated Analyser.

**IDS-iSYS Cuvettes Cube** : IS-CC100, pack of 1 x 960 cuvettes.

**IDS-iSYS System Liquid** : IS-CS100, 5L bottle, ready to use.

**IDS-iSYS Wash Solution** : IS-CW100, 10L bottle, ready to use.

**IDS-iSYS Trigger Set** : IS-CT100, 2 x 250ml vial ready to use.

REF IS-6010 	IDS-iSYS CARTRIDGE CHECK SYSTEM (CCS)	 immundiagnosticsystems
Package Insert	IN VITRO DIAGNOSTIC	

## Assay Procedure

### Preparation and use of CCS on the IDS-iSYS analyser

Proceed according to the instructions of the IDS-iSYS user operational manual.

To have an IDS-iSYS analyser fully operational at the beginning of the day, it is recommended to load the CCS cartridge on board the analyser the evening before, then request the analyser to wake up at a suitable hour in the morning, which will automatically start the analyser performance check using the CCS cartridge.

The reagents provided in the cartridge are ready to use. The analyser automatically performs mixing of the magnetic particles until homogeneity.

The barcode must be read when the cartridge is loaded on the reagent tray. A manual procedure exists to enter the barcode data if the label cannot be read by the analyser barcode reader (see User manual).

Load the cartridge on the reagent tray and wait for 15 minutes before starting the assay.

If there is insufficient reagent present in the cartridge and a second cartridge is required to complete the CCS profiles, load another cartridge of the same lot ensuring the new cartridge is placed on a higher position number than the cartridge in use.

If the cartridge has to be removed from the reagent tray, cap the container with Parafilm and store the cartridge vertically at 2 – 8 °C.

Check the proper identification of the cartridge by the user interface (GUI), record it in the analyser database by reading the mini CD provided.

### Performing the CCS on the IDS-iSYS analyser

All of the data needed to use the CCS cartridge batch can be found on the mini CD.

Start the performance check of the analyzer according to the IDS-iSYS user operational manual.

The measurements requested by the performance check are performed in triplicate for the CCS1, CCS2, CCS3<sup>1</sup>, CCS3<sup>2</sup>, CCS3<sup>3</sup> and CCS3<sup>4</sup> (value in superscript indicates the washer number used in the CCS3 testing).

### CCS frequency:

The analyser performance check is recommended to be run daily, either at the beginning or end of the working day.

A new set of CCS tests is recommended:

- When the QC data does not fall within the requested/defined ranges;
- After instrument maintenance.

### Data processing

The mean, standard deviation, minimum and maximum values, and coefficient of variation are computed from the raw data measured in the various CCS tests. This data can be analysed using the Excel spreadsheet on the mini CD.

## Validation criteria of the CCS

The validation criteria for the CCS are available on the QC certificate on the mini CD-ROM for each lot of CCS

## Procedure in case invalid CCS results

Refer to the user manual of the IDS-iSYS analyser.

IDS S.A.: Rue E. Solvay 101, 4000 Liège, Belgium  
Tel : +32 (0)4 252 26 36, Fax : +32 (0)4 252 51 96  
e-mail : info.be@idsplc.com • www.idsplc.com

**UK** Immunodiagnostic Systems Ltd (IDS Ltd), 10 Didcot Way, Boldon Business Park, Boldon, Tyne & Wear, NE35 9PD  
Tel: +44 (0) 191 519 0660 • Fax: +44 (0) 191 519 0760  
e-mail: info.uk@idsplc.com • www.idsplc.com

**USA** Immunodiagnostic Systems Inc (IDS Inc.), 8425 N. 90th Street, Suite 8,, Scottsdale, AZ 85258  
Tel: 480 278 8333 • Fax: 480 836 7437  
e-mail: info.us@idsplc.com • www.idsplc.com

**Germany** Immunodiagnostic Systems GmbH (IDS GmbH), Mainzer Landstrasse 49, 60329 Frankfurt am Main  
Tel: +49 (0) 69 3085-5025 • Fax: +49 (0) 69 3085-5125  
e-mail: info.de@idsplc.com • www.idsplc.com

**France** Immunodiagnostic Systems SA (IDS SA), 153 avenue d'Italie, 75013 Paris  
Tel: +33 (0)1 40 77 04 50 • Fax: +33 (0)1 40 77 04 55  
e-mail: info.fr@idsplc.com • www.idsplc.com

**Scandinavia** Immunodiagnostic Systems Nordic a/s (IDS Nordic a/s), Marielundvej 30, 2. Sal, 2730 Herlev, Denmark  
Tel: +45 44 84 0091  
email: info.nordic@idsplc.com • www.idsplc.com

REF IS-6010 	SYSTEME DE CONTROLE DE L'IDS-iSYS PAR CARTOUCHE (CCS)	 Immunodiagnostic systems
Notice	DIAGNOSTIC IN VITRO	CE

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Immunodiagnostic Systems S.A., Rue E. Solvay 101, B-4000 Liège, Belgien
**USA** Immunodiagnostic Systems Inc (IDS Inc.), P.O. Box 17063, Fountain Hills, AZ 85269-7063
**Tyskland** Immunodiagnostic Systems GmbH (IDS GmbH), Mainzer Landstrasse 49, 60329 Frankfurt am Main
**Belgien** Immunodiagnostic Systems S.A., Rue E. Solvay 101, 4000 Liège